



Rise & Recline Ltd.

"We adjust to your needs"



FOR EASY LIVING

BESPOKE RISE AND RECLINE CHAIRS

HANDCRAFTED ADJUSTABLE BEDS


QUALITY ACCESSORIES





CONTENTS

P3.	THE PHILOSOPHY
P4.	INTRODUCTION - CHAIRS
P5.	REGAL
P6.	ROYAL
P7.	HAMPTON
P8.	SANDRINGHAM
P9.	T-BACK
P10.	LATERAL
P11.	NOBLE
P12.	BUCKINGHAM
P13.	KENSINGTON
P14.	SETTEES
P15.	BARROWFORD
P16.	MARIANNE
P17.	INTRODUCTION - BEDS
P18.	BALMORAL
P19.	BLENHEIM
P20.	WINDSOR
P21.	RICHMOND
P22.	CROMWELL
P23.	BEDLIFTERS
P24.	HEADBOARDS AND MATTRESSES
P25.	CHAIR AND BED ACCESSORIES
P26.	INTRODUCTION - MEDICAL
P27.	CLARENCE
P28.	FLORENCE
P29.	CAVELL
P30.	VERSATILITY
P31.	REVIVE
P32.	BARIATRIC
P33.	BARROWFORD
P34.	CLARA
P35.	SPECIFICATIONS AND SELECTION TABLE



ATTENTION TO DETAIL, UNSURPASSED QUALITY, SUPERIOR SERVICE.
THE RISE & RECLINE PHILOSOPHY BEGINS HERE...

Rise & Recline Ltd. has been successfully servicing the needs of the mobility industry since 2001 through a combination of superbly designed and crafted chairs, beds and accessories, rapid turnaround of orders and a dealer support service that is second to none.



Rise & Recline Ltd. has always believed in giving its customers the products that best suit their individual needs. This means listening closely to its dealers to identify new market trends and, where applicable, responding quickly to develop new products that can be recommended to their customers with confidence.

Through its state-of-the-art facility in Long Eaton, Nottingham, Rise & Recline Ltd. is demonstrating its commitment to upholding the traditional virtues of British design and craftsmanship. This philosophy runs through every item it produces and is the reason Rise & Recline Ltd. continues to set the highest standards in the industry.



THE CHAIRS

A TRULY BESPOKE RANGE HANDCRAFTED TO MEET EVERY REQUIREMENT

The range of bespoke adjustable chairs and settees defines the market for British manufactured, high quality, versatile seating. Whatever your customer's age, size or ability, Rise & Recline Ltd. is confident you will not find a higher quality chair to match the customer's individual requirements now and in the future.

Please note that when selecting a chair with your customers you should first choose the mechanism, then the style and finally the fabric. The choice of fabric is unlimited to suit any taste or to match existing items of furniture.

The following mechanisms can be incorporated into any style of Rise & Recline Ltd. adjustable chair:

WALL HUGGER

Single motor mechanism that requires only 4-5" of clearance behind to recline.

DUAL MOTOR

Back and footrest are operated independently by separate motors.

2 WAY TILT-IN-SPACE SINGLE

Single motor mechanism that maintains angle of back and seat when reclining and lifting.

2 WAY TILT-IN-SPACE DUAL

A 2 way tilt-in-space with the added benefit of a motorised independent back movement.

TRIOLIFT™

2 way tilt, vertical lift and combination single or dual motor action for ease of transfer and carer contact.

RECLINE ONLY

Single motor wall hugger action mechanism that allows recline only.

MANUAL

Manually operated wall hugger recline mechanism for left or right handed operation.

HEAVY DUTY

Single and dual motor mechanisms available for users up to 40 stone.

BARIATRIC

Three motor assembly to power footrest, back and lift on chairs for customers between 50 and 70 stone.



THE REGAL

SUMPTUOUS STYLE,
TIMELESS BEAUTY



With its small delicate wings and broad, strikingly beautiful button back, there is no chair more inviting than the Regal.

Traditional elegant styling combined with timeless good looks allow the Regal to blend in effortlessly with most home environments.



Removable arm caps made in matching fabric. See accessories on page 25.

Specifications

- 2 pockets
- Duvet seat
- Heavy duty option available
- Handset loop
- Emergency lowering

Accessories - see page 25

Practicality and versatility combined with modern styling are the key features which make the Royal one of the most popular chairs on the market. Large supportive wings ensure the cushions are held in position while the three bolster cushions, each removable with its own zipped back, enable the users to create their own personalised configuration.

THE ROYAL

BESPOKE COMFORT,
FLEXIBLE SUPPORT



Optional removable, swing or lift up arms ensure side transfers are as quick and simple as possible. See accessories on page 25.

Specifications

- 2 pockets
- Duvet seat
- Heavy duty option available
- Handset loop
- Emergency lowering

Accessories - see page 25



THE HAMPTON

TRADITIONAL STYLING,
SIMPLISTIC ELEGANCE



The Hampton boasts all the key features of the Regal but with larger wings to create a back that looks both neat and elegant. Traditional styling and timeless good looks, complete with a beautiful button back, enable the Hampton to complement a wide range of traditional three piece suites.

Specifications

- 2 pockets
- Duvet seat
- Heavy duty option available
- Handset loop
- Emergency lowering

Accessories - see page 25

THE SANDRINGHAM

CLASSIC AND COMFORTABLE WITH A CONTEMPORARY FEEL.



A classic looking removable flat fibre filled back with large wings. The back style on this chair gives the chair a more modern and contemporary feel.

Specifications

- 2 pockets
- Duvet seat
- Heavy duty option available
- Handset loop
- Emergency lowering

Accessories - see page 25



THE T-BACK

VERSATILE DESIGN, SPECIALISED SPINAL SUPPORT



A versatile, three pillowed back with removable cushions. This provides a gap down the centre which is particularly beneficial to customers with spinal problems as pressure is relieved when sitting for long periods.

Specifications

- 2 pockets
- Duvet seat
- Heavy duty option available
- Handset loop
- Emergency lowering

Accessories - see page 25

THE LATERAL

SUPPORTIVE DESIGN,
EASILY ADAPTABLE TO
INDIVIDUAL REQUIREMENTS



The supportive, fibre filled pillow back with removable panels provides superb comfort and the velcro attachments allow for individual adjustments to be made. The filling can be personalised to meet the customer's specific requirements, as the panels have zipped access.

Specifications

- 2 pockets
- Duvet seat
- Heavy duty option available
- Handset loop
- Emergency lowering

Accessories - see page 25

THE NOBLE

CONTEMPORARY AND
FUNCTIONAL

A modern flat foam filled back with no wings.
This firm back style gives this Rise & Recline chair a more contemporary and fresh look.



Specifications

- Duvet seat
- Heavy duty option available
- Handset loop
- Emergency lowering

Accessories - see page 25





THE BUCKINGHAM

MODERN, STYLISH DESIGN,
COMPLEMENTS MOST ROOM SETTINGS



The Buckingham is a modern, stylish rise and recline chair, making it an excellent choice, as it also complements many traditional settees and armchairs. The Buckingham can incorporate all the extra features available with any of the adjustable chairs from Rise & Recline Ltd.



Specifications

- Duvet seat
- Integrated handset
- Heavy duty option available
- Emergency lowering

Accessories - see page 25

THE KENSINGTON

CONTEMPORARY DESIGN,
LUXURIOUS COMFORT



The perfect combination of comfort and style, the Kensington boasts modern and contemporary design with opulent and comfortable cushioning on both the back and the arms. The controls are integrated into the arms of the chair, ensuring that the Kensington complements your standard living room furniture. Tilt-in-Space functionality not available on this style.

Specifications

- 2 pockets
- Integrated handset
- Duvet seat
- Emergency lowering
- Heavy duty option available

Accessories - see page 25





Rise & Recline Ltd. will handcraft adjustable settees to match customers' existing furniture or to complement other chairs.

Handcrafted to complement any of the nine styles of chairs that Rise & Recline Ltd. manufacture, the two and three seater settees provide the perfect solution for those who enjoy sitting together. Settees are available in three versions – static, manual recline or electric recline.

Specifications

- 2 pockets
- Handset loop
- Duvet seat
- Emergency lowering

Accessories - see page 25



The Barrowford offers a variety of styles which can be made with or without side infills and/or wings. A Royal or Lateral Barrowford is also available.

A classic fireside chair, fully upholstered and supplied with the option of wings and infills. Customers can specify wood finish. Also available is a matching two seater settee and static or adjustable footstool.

SETTEES

INDIVIDUAL COMFORT,
TOTAL CONTROL



THE BARROWFORD

CLASSIC LOOK,
BEAUTIFULLY CRAFTED

Specifications

- Removable seat cushion

Accessories - see page 25



THE MARIANNE

FLEXIBLE SPECIFICATION,
ROBUST DESIGN



Specifications

- Up to 30" seat width available
- Removable seat cushion

Accessories - see page 25



The Marianne is available with or without wings.

Complementing the heavy duty range of static and electric armchairs, the Marianne fireside chair is a cleanly styled armchair that can be tailored to the individual needs of each customer. The robust design can handle up to 35 stone with an extra wide seat and is upholstered in the customer's choice of fabrics.



THE BEDS

EVERY BED IS CUSTOM MADE TO ENSURE
ULTIMATE COMFORT

The Rise & Recline Ltd. philosophy of superb quality, craftsmanship and attention to detail is evident again throughout its superb range of luxury beds, headboards and accessories.

Each bed is manufactured to the customer's individual requirements in sizes ranging from 2' 3" to 6' 0" wide and 5' 8" to 7' 0" long (internal). You can also specify bed heights between 16" (foam mattress only) and 30". Standard upholstery is high quality cream damask with a choice of other fabrics available.

All beds can be supplied with a choice of four mattresses – pocket sprung, luxury pocket sprung, Reflex foam and memory foam.

HEAVY DUTY

All Rise & Recline Ltd. standard beds (excluding Cromwell) can be upgraded to heavy duty specification with the addition of a special mechanism, which incorporates two sets of motors, giving a total lift capacity of 35 stone. In this way luxury and comfort are not compromised and the look of the bed remains the same.

THE BALMORAL

ULTIMATE LUXURY,
LAVISH DETAILING



This luxurious bed is designed for those customers who appreciate uncompromising quality. The Balmoral, available in single, double and dual, features a fully padded and upholstered surround with integrated head and footboard.

Specifications

- Single widths: 2' 3", 2' 6", 3' 0"
- Double widths: 3' 6", 4' 0", 4' 6"
- Dual widths: 4' 6", 5' 0", 6' 0"
- 6 button handsets
- 5 part, 4 movement actions
- Up to 18 stone capacity motor as standard
- Furniture castors
- Mattress clip
- Piped finish on damask cover
- Up to 35 stone user upgrade available
- Sectional

Accessories - see page 25



The Relaxor 5 motor massage controller. See accessories on page 25.



THE BLENHEIM

LUXURIOUS COMFORT,
STYLISH DESIGN

The Blenheim exudes luxury and comfort, with a stylish and contemporary feel. Available in single, double and dual, the Blenheim features a fully padded and upholstered surround with integrated headboard for maximum comfort. This modern bed is upholstered in a vast range of fabrics as standard to complement any style of bedroom.



VRT massage
Infrared remote handset containing both massage and bed controls. Night light is available as standard on both single and double zone systems.



Specifications

- Single widths: 2' 3", 2' 6", 3' 0"
- Double widths: 3' 6", 4' 0", 4' 6"
- Dual widths: 4' 6", 5' 0", 6' 0"
- 6 button handsets
- 5 part, 4 movement actions
- Up to 18 stone capacity motor as standard
- Furniture castors
- Mattress clip
- Piped finish on damask cover
- Up to 35 stone user upgrade available
- Sectional

Accessories - see page 25

THE WINDSOR

SUMPTUOUS STYLE,
MAXIMUM VERSATILITY



Designed for less mobile customer, the highly popular Richmond combines luxury and supreme comfort with practicality. Its double arch design permits easy hoist and over-bed table access.

The Richmond comes complete with a padded surround, decorative cording and a stylish footboard.



With its totally padded surround the Windsor offers the perfect divan style option. The versatile design allows for the optional fitting of either one or two drawers in the end or side sections, depending on the width of the bed. This provides valuable additional storage space.

Specifications

- Single widths: 2' 3", 2' 6", 3' 0"
- Double widths: 3' 6", 4' 0", 4' 6"
- Dual widths: 4' 6", 5' 0", 6' 0"
- 6 button handsets
- 5 part, 4 movement actions
- Up to 18 stone capacity motor as standard
- Furniture castors
- Mattress clip
- Piped finish on damask cover
- Up to 35 stone user upgrade available
- Sectional

Accessories - see page 25

The Windsor side drawer option. See accessories on page 25.



THE RICHMOND

ULTIMATE COMFORT,
SPECIALIST DESIGN

Grab rails are available for added support when users are entering or leaving the bed. See accessories on page 25.



Specifications

- Single widths: 2' 3", 2' 6", 3' 0"
- Double widths: 3' 6", 4' 0", 4' 6"
- Dual widths: 4' 6", 5' 0", 6' 0"
- 6 button handsets
- 5 part, 4 movement actions
- Up to 18 stone capacity motor as standard
- Furniture castors
- Mattress clip
- Piped finish on damask cover
- Up to 35 stone user upgrade available
- Sectional

Accessories - see page 25



A traditionally styled single bed designed with practicality and affordability in mind. The Cromwell is supported by four polished wooden skittle legs and, with 8" clearance, is ideally suited for hoists and over-bed table access. Birch hardwood support rails and corner gussets ensure a strong but simple design, complemented by high quality upholstery and a choice of mattresses.

Specifications	
• Single widths: 2' 3", 2' 6", 3' 0"	
• 6 button handset	• Mattress clip
• 5 part, 4 movement action	• Up to 18 stone capacity motor as standard
Accessories - see page 25	



THE CROMWELL

OUTSTANDING VALUE,
UNCOMPROMISING
FLEXIBILITY AND QUALITY



BED LIFTERS

ASSISTANCE MAXIMISED,
AESTHETICS MAINTAINED



HI-LO PROFILER
This bed lifter is incorporated into the bed surround and lifts the mattress and profiling mechanism only to a maximum height of 38" (including mattress). Available in 2' 6" and 3' 0" widths, it is perfect for single beds and for dual beds where couples wish to retain their closeness at night without losing the practicality of a care bed.

VERTICAL BED LIFTER
This high quality mechanism attaches to the bed surround and lifts the entire assembly vertically to a maximum height of 41" (including mattress). Four lockable heavy duty castors enable ease of movement and ensure that the 5" gap for hoist access is maintained. Gas springs assist in the smooth operation of the mechanism during raising and lowering. OKIN actuators, along with an in-line transformer fitted with a mains isolating plug, ensure that power to the lifting mechanism is only supplied when the handset is activated. An emergency recovery pack is also fitted. Available in widths from 3' 0" to 4' 6".



HEADBOARDS

Designed and handcrafted in the Rise & Recline Ltd. workshops, the luxuriously padded and beautifully detailed headboards are the ideal complement to the Cromwell, Richmond and Windsor beds.



Shell – Elegantly sculpted shell design for extra flair
D/D – Classic design with single plain border

Fraser – Contemporary borderless design with piped edging

MATTRESSES


Pocket sprung
A traditional 8” sprung mattress manufactured with comfort and support in mind. Individually pocketed springs are designed to follow the contours of an adjustable bed. Upholstered in a matching cream damask fabric with brass ventilators and side handles for ease of lifting and turning.

Luxury pocket sprung
Similar to the pocket sprung but with extra padding and a greater thickness (9”) for additional comfort and support.

Reflex
Made from high quality foam and finished in a superior quality COOLMAX® two way stretching fabric, the Reflex 8” mattress has an outstanding look and feel. The foam is available in various densities for different user weights and ensures a comfortable and relaxing sleep.

Memory foam
The 8” memory foam mattress combines both Reflex foam and memory foam, so providing pressure relief and enabling the correct posture to be maintained during sleep. Finished in a superior quality COOLMAX® two way stretch fabric, the memory foam mattress cradles the body for the ultimate in comfort and relaxation.

COOLMAX®
The standard cover for the Reflex or memory foam which draws heat out of the mattress to avoid overheating.



- FOAM COLOUR KEY
- 0-12st – pink 400P
 - 12-15st – grey 400H
 - 15-18st – white 400X
 - 18st+ – white 650/F

MASSAGE OPTIONS
Both the Reflex and memory foam mattresses are available with an inbuilt Relaxor 5 or Relaxor 10 motor massage system – see accessories on page 25.

Accessories - see page 25

...FOR CHAIRS

Swing Arm

Removable Arm

Lift Up Arm

Side transfer arms - swing, removable or lift
This option enables side transfers out of the chair with the minimum of difficulty. Left hand, right hand or both.

Removable arm caps
Made from matching fabric, these removable covers help protect the chair arms from wear and tear.

Arm knuckles
High quality polished wood knuckles for extra leverage and convenience.

Head pillow
For extra comfort and support. Made in matching fabric.

Pressure seating
Refer to separate Specialist Seating Guide.

Head cover
Protective fabric cover suitable for all chairs.

Fabric protection
All chairs can be treated with McKLords protection to help maintain the fabric's quality and appearance.

Heat system
Optional heat pads in the seat or back or both. Simple two button hand control.

6/10 motor massage
The Relaxor massage system has either six or ten motors built into the chair and comes with a heat pad for the lumbar region. A hand controller allows the user to adjust the speed and intensity of each zone.

Removable seat
Versatile cushion cut out that allows pressure care to be changed easily.

Duvet style removable seat
Versatile cushion cut out that allows pressure care to be changed easily housed underneath the duvet.

ACCESSORIES

...FOR BEDS

5/10 motor massage
The Relaxor 5 massage system has five motors integrated equally within the mattress (foam only). The Relaxor 10 massage system has ten motors integrated equally in five zones within the mattress to give a more powerful massage (foam only). Normally used in 4' 6", 5' 0" and 6' 0" double mattresses.

Cot sides
Three or four bar chrome plated steel side rails which contain the occupier safely within the bed area.

Divan drawers
Available for the Balmoral, Blenheim and Windsor beds. Can be situated at the end or on the side.

VRT massage
Infrared remote handset containing both massage and bed controls. Night light is available as standard on both single and double zone systems.

Grab rails
Solid steel construction to help occupant entering or leaving the bed.

Length alteration
All beds have the option of being shortened or lengthened. Internal sizes available are 5' 8", 6' 0", 6' 3", 6' 9" and 7' 0". Standard length = 6' 6".

Height alteration
All beds have the option of the total height (floor to top of mattress) being raised to 30" or lowered to 16" depending on the mattress. Reflex/memory foam from 16-17", pocket sprung/luxury pocket sprung from 18-19".

Slide back
Slides the whole action backwards like a wall hugger mechanism. Available in 2' 6" and 3' 0".

Mattress covers
Available as fitted sheets or permanent waterproof covers. Sizes from 2' 6" to 6' 0".

Link twin option
Links together two single beds to form a dual.

- NHS • Local Authority • Nursing Homes •

INSPIRED BY CLINICIANS AND CREATED BY US

Rise & Recline Ltd. is a strong provider to the NHS, Local Authority and nursing homes, who have earned the respect of key clinicians within the industry. By working in partnership, they have harmonised main stream products into the industry to meet the patient's varying needs.



MEDICAL

Rise & Recline Ltd. is confident that, if you have the desire for inspiration, innovation and change they have the skills and expertise to deliver and meet your expectations.

Rise & Recline Ltd. already offer a wide range of innovative seating products, created by need, without compromise on patient safety, ensuring a free from harm environment. They understand that different professionals require diverse expectations from any postural supported seat.

Rise & Recline Ltd. work alongside infection control, tissue viability and moving and handling departments to ensure the ergonomics and pressure relief characteristics are met and maintained.

Rise & Recline Ltd. is one of the finest manufacturers of general, bariatric and paediatric seating in the UK paired with a customer service that is second to none.



Rise & Recline Ltd.
"We adjust to your needs"



SAFE WORKING LOAD
• 160kg / 25st

THE CLARENCE

VERSATILE, COMFORTABLE CHAIR, TILT-IN-SPACE ACTION, ELECTRIC OR MANUAL

Specifications

- Velcro lateral back
- Available as manual
- Removable seat cushion
- Hoist accessible

VARIED SEAT HEIGHTS, WIDTHS AND DEPTHS

- Ergonomically designed to accommodate most body shapes

PRESSURE REDUCING CUSHIONS

- Allows an even distribution of pressure
- May reduce the risk of pressure ulcers
- Maintains skin integrity
- Increases patient comfort

INDEPENDENT POWERED LEG LIFT

- May increase circulation
- Supports the patient when raising the lower limbs



POWERED TILT-IN-SPACE ACTION

- Limited shear and friction
- Promotes micro movements and may reduce the risk of pressure ulcers
- Promotes weight distribution
- Maximises patient comfort and support

CASTORS

- Low friction castors
- 3 Locking castors (1 Directional) and 1 standard

HOIST FRIENDLY

- Aids in patient transfers
- Allows precise patient positioning

DURABLE MATERIAL

- Waterproof, which inhibits the growth of bacteria
- Reduces tissue loading, shear and friction
- Extremely durable and promotes longevity





FLORENCE

ERGONOMICALLY
DESIGNED WITH
PRACTICALITY IN MIND

Specifications

- 4 x large low friction, multi directional braked castors
- SWL 25 stone / 160kg
- Multi positional, manual self-locking leg support
- 5 year guarantee
- High tech durable material
- Hygiene access
- Ergonomically positioned hand bar
- Pressure reducing cushions
- Elegant simple design



BENEFITS

- Can be transferred around the home, across thresholds and very easy to manoeuvre
- Promotes micro movements, maximises patient comfort and independence, heightens patient support, can improve circulation
- Risk adverse for 5 years
- Waterproof, Anti-bacterial / effective against C.Difficile spores, bleach cleanable, Crib 5 flame retardant
- Allows any excess fluids to escape and easy for maintenance
- Enables the carer to manoeuvre the chair with ease whilst maintaining a neutral position
- Allows even distribution of pressure, assists in the prevention of pressure ulcers, maintains skin integrity, increases patient comfort
- Will cater for 90% of your patients weights and sizes
- Fixed back seat angle
- Optimum postural seating position



BENEFITS

- Enables smooth transfers around the room and across thresholds
- Ideal for patient room to room transportation and very easy to manoeuvre
- Risk adverse for 5 years
- Waterproof, Anti-bacterial / effective against C.Difficile spores, bleach cleanable, Crib 5 flame retardant
- Allows any excess fluids to escape and easy for maintenance
- Enables the carer to manoeuvre the chair with ease whilst maintaining a neutral position
- Allows even distribution of pressure, assists in the prevention of pressure ulcers, maintains skin integrity, increases patient comfort
- Will cater for 90% of your patients weights and sizes
- Manual 2-Way Tilt-in-Space recline function
- Optimum postural seating position



Specifications

- 4 large low friction, multi directional braked castors
- SWL 25 stone / 160kg
- Available with Royal waterfall or Lateral cushions
- 5 year guarantee
- High tech durable material
- Hygiene access
- Ergonomically positioned hand bar
- Pressure reducing cushions
- Elegant simple design
- Multi positional, manual self-locking leg support, with robust bumper protection



CAVELL

DURABLE, PRACTICAL
AND SIMPLE

THE VERSATILITY

INGENIOUS AND SIMPLE
ADJUSTABLE, MODULAR
DESIGNED CHAIR



BENEFITS

- Varied seat heights, widths and depths in one chair
- Freedom to change pressure reducing cushions easily and quickly
- Varied Mechanism options available
- Tilt-in-Space, Tilt-in-Space Dual, Dual.
- Modular “tool-less” design to enable efficient and effective repair, reducing cost and increasing recyclability
- Seating/postural solutions and vast accessories also available with this concept
- Various Upholstery finishes and colours to suit a multitude of environments



**Removable Wings For
Optimum Recyclability**



Adjustable Postural Support



Adjustable Pillow Backrest
(various styles available)



Removable Cushioning



**Integrated Width Adjustment
System**
(patented design)



Fully Removable Duvet



Height Adjustment System



THE REVIVE

MODULAR DESIGNED
CHAIR, PERFECT FOR LOAN
STORE SUPPLY



BENEFITS

- Freedom to change pressure reducing cushions easily and quickly
- Varied mechanism options available – 2-Way Tilt-in-Space Single or Dual
- Modular “tool-less” design to enable efficient and effective repair, reducing cost and increasing recyclability
- Seating/postural solutions and vast accessories also available with this concept
- Various Upholstery finishes and colours to suit a multitude of environments



**Removable Wings For
Optimum Recyclability**



Fully Removable Duvet



Fully Removable Headset

BARIATRIC

SPECIALIST DESIGN,
POWERFUL MOVEMENT



VARIED SEAT HEIGHTS, WIDTHS AND DEPTHS

- Ergonomically designed to accommodate most body shapes
- Allows maximum patient support, safety and comfort

PRESSURE REDUCING CUSHIONS

- Allows an even distribution of pressure
- May reduce the risk of pressure ulcers
- Maintains skin integrity
- Increases patient comfort

POWERED LEG LIFT

- Increases circulation
- May reduce under-shearing
- Supports the patient when raising the lower limbs

POWERED BACK REST

- Maximises patient comfort and support
- Promotes micro movements and reduces the risk of pressure ulcers
- Promotes weight distribution

HOIST FRIENDLY (Overheads/Gantry)

- Aids in patient transfers
- Allows precise patient positioning

DURABLE MATERIAL

- Waterproof, which inhibits the growth of bacteria
- Reduces tissue loading, shear and friction
- Extremely durable and promotes longevity



SAFE WORKING LOAD

- 260kg / 40st to 450kg / 70st
- 22st Leg Lift
- Perfect to suit most bariatric/plus size groups



Specifications

- | | |
|--------------------|----------------------|
| • 2 pockets | • Handset loop |
| • Duvet seat | • Emergency lowering |
| • Pressure seating | • Heavy duty foam |

BARROWFORD

CLASSIC LOOK,
BEAUTIFULLY CRAFTED



Specifications

- SWL 25 stone / 160kg
- Crib 5
- High tech durable material
- Pressure reducing cushions
- Hygiene access
- Ergonomically positioned hand bar and housekeeping wheels available on request
- Air cushions
- Elegant simple design
- Aesthetically pleasing



BENEFITS

- Ideal for all bedside and communal locations in the NHS, nursing homes and end user markets
- Waterproof, Anti-bacterial / effective against C.Difficile spores, bleach cleanable
- Inter changeable seat cushion for medium to high risk patients, without changing the seat height
- Enables the carer to manoeuvre the chair with ease around the room and wards
- Will cater for an alternating air cushion and pump without compromising the ergonomics of the chair
- Fixed back seat angle
- Optimum postural seating position
- Can be made in various colours to blend into most environments
- Will cater for 99% of patients, with available upgrade to 35 stone



THE CLARA

A HIGH QUALITY, DURABLE
AND FEATURE RICH
MEDICAL BED



BENEFITS

- 4 Section Electric Profiling
- IEC 60601 – 2-52 Compliant
- Trendelenburg tilt to -12°, Reverse Trendelenburg tilt to +12°
- Perceptive Profiling System (PPS) minimises Friction & Shear
- Soft Synthetic Castors (100mm)
- Integral Bed Extension of 17.5cm
- Lock-out feature handset
- GoLow GoSlow safety feature
- Vascular and Fowler leg rest positioning
- Curved head and foot board design provides patients with 180° visual field
- Enforced PU deck slats offer improved infection control and reduce lifetime maintenance costs
- Anti-Bacterial powder coated frame

Specifications

- 4 Section Electric Profiling up to 250kgs
- Trendelenburg tilt to -12°, Reverse Trendelenburg tilt to +12°
- GoLow GoSlow Safety Feature
- Infection Control
- 250kg weight capacity
- Vascular and Fowler leg rest positioning
- Standard Bed Height Range 23.5cm to 69cm
- Low Bed Height Range 13.5cm to 59cm - ideal for use with a cushioned safety mat



CHAIR SELECTION

Mechanism	Weight Capacity	Seat Height	Seat Depth	Seat Width	Trio
Static	SWL 20st	15" – 21"	17" – 21"	18" – 28"	
Manual Recline	SWL 18st	18" – 21"	17" – 21"	18" – 22"	
Recline Only	SWL 18st	18" – 21"	17" – 21"	18" – 22"	
Single Motor Wall-hugger	SWL 18st	18" – 21"	17" – 21"	18" – 22"	
Dual Motor (2)	SWL 18st	18" – 21"	17" – 21"	18" – 22"	
Low Profile 3-Way Single	[SWL 20st]	[15" – 17"]	[17" – 21"]	[18", 20" & 21"]	
Low Profile Dual (2)	SWL 20st	15" – 17"	17" – 21"	18" – 20" [21" – 22"]	
	[SWL 35st]	[15" – 17"]	[17" – 21"]	[22" – 28"]	
Heavy Duty 3-Way Single	[SWL 25st]	[18" – 21"]	[17" – 21"]	[24", 26" & 28"]	
	[SWL 35st]	[18" – 21"]	[17" – 21"]	[24", 26" & 28"]	
Heavy Duty Dual	[SWL 25st]	[18" – 21"]	[17" – 21"]	[22" – 28"]	
	[SWL 35st]	[18" – 21"]	[17" – 21"]	[22" – 28"]	
2-Way Tilt-in Space Single	SWL 20st	17" – 21"	17" – 21"	18" & 20"	
	SWL 20st	15" – 17"	17" – 21"	18", 20" & 22" [24"]	
	SWL 25st	17" – 21"	17" – 21"	18", 20" & 22" [24", 26" & 28"]	Yes
	[SWL 35st]	[17" – 21"]	[17" – 21"]	[22", 24", 26" & 28"]	
2-Way Tilt-in Space Dual	SWL 20st	17" – 21"	17" – 21"	18" & 20"	
	SWL 20st	15" – 17"	17" – 21"	18", 20" & 22" [24"]	
	SWL 25st	17" – 21"	17" – 21"	18", 20" & 22" [24", 26" & 28"]*	Yes

BLACK - 7 - 10 day turnaround
[RED] - Available with longer lead time and price supplement

BED SELECTION

0–18st (0–120kg)	
Cromwell	Single: 2'3"–3'0"
0–35st (0–220kg)	
Balmoral	Single: 2'3"–3'0" Double: 3'6"–4'6" Dual: 4'6"–6'0"
Blenheim	
Windsor	
Richmond	
35–50st (220–320kg)	
Bariatric (hospital style bed)	Single: 3'6"–5'0"

MEDICAL SELECTION

Bariatric		Clarence	
3 motor (40 & 50st) (dual operation)	18–21" seat height 18–21" seat depth 22–30" seat width	Manual 2-way tilt-in-space	11–20" seat height (adjustable footboard) 16–23" seat depth 16–21" seat width
3 motor (50 & 70st) (dual operation)	18–22" seat height 19–21" seat depth 30–34" seat width	Electric 2-way tilt-in-space (manual back)	
Florence		Cavell	
Tilt position transfer chair with ratchet footboard	20" seat height 20" seat depth 20" seat width	Manual Tilt-in-Space Recline with Ratchet footboard	20" seat height 18" seat depth (to end of KB = 21") 20" seat width
Barrowford		Clara	
Fireside static chair	16½–20" seat height 16–20" seat depth 18–20" seat width	Standard bed	80cm bed width 198cm bed length 196cm mattress
		Low bed	

HD = Heavy Duty

Variations may occur to the colour of fabrics and products featured in this brochure as a result of the printing process.

Due to a policy of continual improvement, Rise & Recline Ltd. reserve the right to change specifications without prior notice.



THE ORIGINAL MANUFACTURER OF RISE AND RECLINE
BESPOKE CHAIRS AND ADJUSTABLE BEDS



Rise & Recline Ltd.
"We adjust to your needs"

t: 0115 913 3572 f: 0115 913 3570 e: info@riseandrecline.co.uk www.riseandrecline.co.uk
Units 5 & 6, Palmer Business Court, Manor House Road, Long Eaton, Nottinghamshire NG10 1LZ

Registered Office: Units 5 & 6, Palmer Business Court, Manor House Road, Long Eaton, Nottinghamshire NG10 1LZ
Company Registered in England and Wales No. 419 4331 . VAT Registered No. 756 0095 30